

## Product News January 2009

### About MEDACX

MEDACX is an established supplier to the UK's First Responders, Paramedics, Ambulance Services, Doctors, Nurses, Hospitals, A&E, Armed Forces Medical and General workplace medical organisations; providing essential portable monitors, sensors and diagnostic devices designed for use within the emergency services environments.



### **EMMA™** Emergency Capnometer

*A must have for every emergency grab bag!*

### Capnography

*EMMA™ is a truly portable mainstream capnometer, using infrared technology providing rapid, accurate and reliable, End tidal CO<sub>2</sub> and Respiratory Rate measurements.*



*EMMA™ connects directly to LMA's, laryngeal masks, endotracheal tubes, or facial masks using the unique disposable airway adapter. Designed for emergency monitoring for use on non-intubated patients, during CPR and for intubated patients, where it helps ensure correct airway placement providing indication of efficacy of ventilation.*



*All products are application proven and come from leading European medical device companies such as; Bluepoint, Defibtech, HF, Phasein and EnviteC-Honeywell. Products meet the European medical device directive, are manufactured and quality assured under ISO 13485, carry CE marks, have full manufacturer's warranties and are fully supported by MEDACX.*

*MEDACX has a complete range of products and those not included in this review can be found on our website, or for more information please contact one of our support staff.*

## Product News *January 2009*



## Pulse Oximeters

### **OxyTrue® SC** **with full 24 month warranty!**

Fingertip pulse oximeter that's **designed to withstand the harshest emergency medical environments**. It's totally encapsulated in a bio-compatible Silicone rubber (non-latex) housing providing a soft but secure fit. Dust and water resistant. Easy to clean ensuring improved patient to patient hygiene. With full colour adjustable orientation displays. Full compliance to ISO 9919:2005



## Pulse Oximetry Sensors

### **SoftCap®**

Rubberboot 3rd generation soft sensors. Silicone (non-latex) rubber flexible construction, bio compatible, with improved comfort and a secure patient fit. Easy to clean, ensures utmost patient to patient hygiene. Available to fit most original manufacturer's equipment (OEM).



### **Bluepoint®**

Disposable SpO2 sensors: Paediatric, Neonates, Infants and Adults. Available to fit most original manufacturer's equipment (OEM).



**Contact us and let us provide a competitive quote for ALL your Pulse Oximetry devices and SpO2 sensor requirements. Proven quality, 'single source supply'.**

## Product News *January 2009*



### **OxiPen®**

*The number ONE choice. A truly portable, accurate and reliable pulse oximeter for spot check monitoring of functional arterial oxygen saturation (SpO<sub>2</sub>) and pulse rate (BPM). Combined with SoftTip sensor it is easy to clean and provides a secure fit with improved patient comfort.*

*Optional additional sensors: infants, paediatric or neonates, protects your investment supporting all patient sizes. Full compliance to ISO 9919:2005*



## Pulse Oximeters

### **OxyTrue® A**

*Handheld pulse oximeter for continuous or spot check monitoring of functional arterial oxygen saturation (SpO<sub>2</sub>) and pulse rate (BPM) of adults, infants, paediatric and newborn patients within hospital, clinical & medical practices, rescue services & emergency patient transportation.*

*The unique Sapphire® technology ensures accurate measurements and coupled with the SoftCap® sensor it delivers improved patient comfort. With silicone rubber cover, colour displays, audible and visual alarm functions, short and long term trend capabilities, data transfer via USB cable and PC software. Optional accessories: carry case, mounting kits and other sensors.*

*Full compliance to ISO 9919:2005*

## Pulse Oximeters

### **OxyTrue® FC**

*Fingertip pulse oximeter for spot check monitoring of functional arterial oxygen saturation (SpO<sub>2</sub>) and pulse rate (BPM).*

*Ideal for spot checks in hospital, clinical practices and homecare.*



## Product News *January 2009*

### **Advance-PAK**

*Defibrillator - Replacement Pads*



#### **Quality Replacement Pads from MEDACX**

*Advance-PAK pads: 2 styles Adult and Infant; are designed to work with the majority of manufacturers defibrillators and individual model types. Individually packaged with leads outside. Made in Germany, CE marked, Non Latex.*

**5050 Advance-PAK**  
*(Schiller/Bruker/GE/Marquette)*

**4040 Advance-PAK**  
*(Welch Allyn)*

**3030 Advance-PAK**  
*(ZOLL)*

**2020 Advance-PAK**  
*(Philips/Laerdal/HP/Agilent)*

**1010 Advance-PAK**  
*(Medtronic Physio-Control)*

*Contact us and let us provide a competitive quote for ALL your Defibrillator Pad requirements. Proven quality, 'single source' supply.*

### **BlueTemp®**

*Temperature Sensors*



#### **Rectal - Oesophageal - Skin Contact**

*For reliable and precise monitoring of body surface and core temperature in adults, infants and neonates. Very Accurate  $\pm 0.10^{\circ}\text{C}$  Measuring Range  $0-50^{\circ}\text{C}$  Autoclavable at  $134^{\circ}\text{C}$ , Latex free and designed to fit most manufacturers original equipment.*

### **LifeLine® Automatic AED**

*Contact us for more details*



*MEDACX has a complete range of products and those not included in this review can be found on our website, or for more information please contact one of our support staff.*