


SpiroTrue® A and SpiroTrue® APC

Hot-wire anemometric flow sensors for ICU ventilators and anaesthesia workstations



» PRODUCTS that care
for you «

*Flow sensors based on hot-wire anemometric principle for
ICU ventilators and anaesthesia workstations*



Specification

- Designed in accordance with 93/42/EEC (Annex I and Annex II)
- Product classification in accordance with Annex IX of 93/42/EEC - class IIa

SpiroTrue® A:

- Compatible to the Dräger Spirolog flow sensor REF 8403735

SpiroTrue® APC:

- Compatible to the Dräger Spirolog REF 8403735, SpiroLife REF MK01900

Applications

- Devices without RFID functionality e.g. Evita family, Savina family, Fabius family, Primus/Apollo, Julian, Cato, Cicero

Note

- SpiroTrue® A/APC flow sensors can be disinfected
- SpiroTrue® APC flow sensor can be autoclaved (134 °C) ≤ 100 cycles

 Made in Germany

Ordering information

SpiroTrue® A:

Order number: 1030132006

Delivery package: SpiroTrue® A flow sensor,
5 pieces per box, individually packaged

SpiroTrue® APC:

Order number: 3030131003

Delivery package: SpiroTrue® APC flow sensor, 1pc/box

For further information visit: www.medacx.co.uk

bluepoint® MEDICAL United Kingdom partner

MEDACX Limited
Units 1-3 Alexander House, 60-62 Station Road
Hayling Island, Hampshire PO11 0EL

Tel: 02392 469737

E-mail: info@medacx.co.uk